



STERILIZATION CONTROL: “BOWIE & DICK PACK”

- Air leaks inside the autoclave represent a serious threat for the validation of the sterilization process. However, this problem can be solved through the daily use of the Bowie & Dick test pack, as provided by international standards in force.
- Disposable pack for the daily Bowie & Dick test is suitable for the control of steam penetration and air removal from the autoclaves during steam sterilization.
- Unlike traditional devices, the use of Bowie & Dick test pack offers consistent advantages that ensure the validation of each test and the immediate detection of air removal during the sterilization.
- Bowie & Dick test pack is a safe and reliable product.
- It can immediately detect any problem occurring during the Bowie & Dick cycle.
- Test sheet will remain flat after the sterilization cycle.
- External chemical indicator permits the prompt identification of processed

REF. EC.STEAM 98



REF. EC.STEAM 98 S



packs.

- The materials are full recyclable and the process indicator is non-toxic.
- Full conformity to standards in force.
- Easy-to-read test sheet clearly identifies failures and can be filed for record keeping.
- 20 pcs/pack.

BOWIE & DICK TEST

REF.	EXTERNAL WRAPPING	REF.	PACK
EC.STEAM 98	cardboard with easy-peel	200 pcs	20 pcs/box
EC.STEAM 98 S	cardboard with easy-peel	160 pcs	20 pcs/box
EC.STEAM 98.1-L	medical grade paper	200 pcs	20 pcs/box
EC.STEAM 98.1/E	creped paper	180 pcs	20 pcs/box
EC.STEAM 98.1/R	use with basket	100 pcs	25 pcs/pack

REF. EC.STEAM 98.1-L



REF. EC.STEAM 98.1/E



STERILIZATION CONTROL: “EC.STEAM 98.1-R” REUSABLE BOWIE & DICK TEST

- Reusable test pack for the control of steam penetration and air removal inside the autoclaves during steam sterilization
- Safe and reliable product in compliance with standards in force
- Each package contains a user-friendly test sheet for data storage
- The kit can be utilized 25 times before disposal
- The kit is composed of two sets of reusable porous sheets, a process indicator for steam monitoring and relative test storage card usable up to 25 times.
- Reusable stainless steel basket with clamp.



REUSABLE BOWIE & DICK TEST